# **CrossMap Documentation**

Release 0.1.3

**Liguo Wang** 

November 20, 2013

# CONTENTS

1	why Crossmap:	3
2	How CrossMap works? 2.1 Algorithm	<b>5</b> 5 5
3	News	7
4	Download	9
5	Installation	11
6	6.2 User Input file	13 13 14 14
7	7.2 Convert BAM/SAM format files	15 16 17 20 21 22
8	Compare to UCSC liftover tool	25
9	LICENSE	27
10	Contact	29



- CrossMap is a program for convenient conversion of genome coordinates (or annotation files) between *dif- ferent assemblies* (such as Human hg18 (NCBI36)  $\Leftrightarrow$  hg19 (GRCh37), Mouse mm9 (MGSCv37)  $\Leftrightarrow$  mm10 (GRCm38)).
- It supports most commonly used file formats including SAM/BAM, Wiggle/BigWig, BED, GFF/GTF, VCF.
- CrossMap is designed to liftover genome coordinates between assemblies. It's *not* a program for aligning sequences to reference genome.
- We do not recommend using CrossMap to convert genome coordinates between species.

CONTENTS 1

2 CONTENTS

**CHAPTER** 

ONE

### WHY CROSSMAP?


Full genome sequencing, especially mammalian (eg. human) genomes, requires extensive, continuous efforts. Therefore reference genome assemblies are subject to change and refinement from time to time. Generally, researchers need to convert results that have been analyzed according to old assemblies to newer versions or *vice versa*, to facilitate meta-analysis, direct comparison as well as data integration and visualization.

Several useful conversion tools have been developed:

- UCSC liftover tool only supports BED input.
- NCBI remap support BED, GFF, GTF, VCF, etc
- Galaxy (Based on UCSC liftover tool) supports BED, GFF, GTF input.
- Ensembl assembly converter supports BED, GFF, GTF, PSL input, but output is GFF only.
- pyliftover "only does conversion of point coordinates, that is, unlike liftOver, it does not convert ranges, nor does it provide any special facilities to work with BED files".

But none have the functionality to convert files in BAM/SAM or BigWig format. This is a significant gap in computational genomics tools, since these formats are the ones most widely used for representing high-throughput sequencing data such as RNA-seq, ChIP-seq, DNA-seq, etc.

#### **HOW CROSSMAP WORKS?**



# 2.1 Algorithm

CrossMap first determines the correspondence between genome assemblies from UCSC chain file (chain file describes the pair-wise alignments between two genomes). Genome intervals will be stored in interval tree data structure, which allows one to efficiently find all intervals that overlap with any given interval or point. Then CrossMap remaps each entry in BAM/SAM, BED, GFF/GTF, VCF file to the target assembly by querying the interval tree. Exon/intron structure in BED file; spliced alignments, paired alignments, insert size, header section, SAM flags in BAM/SAM file; reference alleles, indels in VCF file will be processed properly.

For Wiggle/BigWig format files, line-by-line computation will be very slow. To increase speed, CrossMap groups consecutive coordinates with the same coverage score into bins (i.e. genomic regions), then remaps those regions one-by-one to the target assembly by querying the interval tree. In other words, Wiggle/BigWig files will be converted into bedGraph format internally, which will be converted into BigWig format (if UCSC's 'wigToBigWig' executable exists and is callable).

# 2.2 Time complexity

Assume there are N lines in the chain file. CrossMap loads the chain file first and process the query file line by line. Thus the space complexity is O(N). For each query region (s,t), it takes  $O(\log N)$  time to locate which chain(s) overlap with s and t. Then it takes  $O(\log N)$  time to search the sorted ungapped alignments in this chain that overlap with s and t and calculate the converted values for s and t in the target assembly. So in total it takes  $O(\log N)$  time to convert one query. The time complexity is  $O(\log N*M)$  to convert M queries.

In practical, the time CrossMap takes increases linearly to the size of input file.

# **CHAPTER**

# **THREE**

# **NEWS**

• 10/23/13 4:16 PM: CrossMap (0.1.3) was released

8 Chapter 3. News

# **CHAPTER**

# **FOUR**

# **DOWNLOAD**

- CrossMap source code
- Test datsets

#### **CHAPTER**

#### **FIVE**

### **INSTALLATION**

#### Prerequisite:

- 1. gcc
- 2. python2.7.\*
- 3. numpy
- 4. cython

#### Download CrossMap program from here:

- \$ tar zxf CrossMap-VERSION.tar.gz
- \$ cd CrossMap-VERSION
- # install CrossMap to default location. In Linux/Unix, this location is like:
- # /home/user/lib/python2.7/site-packages/
- \$ python setup.py install
- # or you can install CrossMap to a specified location:
- \$ python setup.py install --root=/home/user/CrossMap
- # setup PYTHONPATH. Skip this step if CrossMap was installed to default location.
- \$ export PYTHONPATH=/home/user/CrossMap/usr/local/lib/python2.7/site-packages:\$PYTHONPATH.
- # Skip this step if CrossMap was installed to default location.
- \$ export PATH=/home/user/CrossMap/usr/local/bin:\$PATH

#### NOTE:

- 1. Due to intensive computation, CrossMap is designed to run on Linux/Unix and Mac OS. Some modules may not work properly on Windows.
- 2. Mac users need to download and install Xcode command line tools.

### **INPUT AND OUTPUT**

CrossMap basically needs 2 input files. chain format file describing genom-wide pairwise alignments between assemblies and the file containing genome coordinates that you want to convert to different assembly. If input file is in VCF format, a reference genome sequence file(in FASTA format) is needed.

#### 6.1 Chain file

#### Example of chain file:

```
chain 4900 chrY 58368225 + 25985403 25985638 chr5 151006098 - 43257292 43257528 1
                  0
         1
 10
         0
                  5
 61
         4
                  0
         0
 16
                  4
         3
 42
 16
 14
         1
                  0
 3
                  0
 48
 chain 4900 chrY 58368225 + 25985406 25985566 chr5 151006098 - 43549808 43549970 2
         0
 60
         4
                  0
 10
         0
                  4
 70
```

UCSC prebuilt most commonly used chain files:

- Human (Homo sapiens)
- hg19ToHg18.over.chain.gz (Chain file needed to convert hg19 to hg18)
- hg19ToHg17.over.chain.gz (Chain file needed to convert hg19 to hg17)
- hg18ToHg19.over.chain.gz (Chain file needed to convert hg18 to hg19)
- hg18ToHg17.over.chain.gz (Chain file needed to convert hg18 to hg17)
- hg17ToHg19.over.chain.gz (Chain file needed to convert hg17 to hg19)
- hg17ToHg18.over.chain.gz (Chain file needed to convert hg17 to hg18)
- Mouse (Mus musculus)
- mm10ToMm9.over.chain.gz (Chain file needed to convert mm10 to mm9)
- mm9ToMm10.over.chain.gz (Chain file needed to convert mm9 to mm10)

• mm9ToMm8.over.chain.gz (Chain file needed to convert mm9 to mm8)

Chain file of other species can be downloaded from http://hgdownload.soe.ucsc.edu/downloads.html

# 6.2 User Input file

- 1. BAM or SAM format.
- 2. BED or BED-like format. BED file must has at least 3 columns ('chrom', 'start', 'end').
- 3. Wiggle format. Both "variableStep" and "fixedStep" wiggle line are supported.
- 4. BigWig format.
- 5. GFF or GTF format.
- 6. VCF format.

# 6.3 Output file

Format of Output files depends on the input format

Input_format	Output_format
BED	BED (Genome coordinates will be updated to the target assembly)
BAM	BAM (Genome coordinates, header section, all SAM flags, insert size will be updated accordingly)
SAM	SAM (Genome coordinates, header section, all SAM flags, insert size will be updated accordingly)
Wiggle	bedGraph (if wigToBigWig executable does not exist)
Wiggle	BigWig (if wigToBigWig executable exists)
BigWig	bedGraph (if wigToBigWig executable does not exist)
BigWig	BigWig (if wigToBigWig executable exists)
GFF	GFF (Genome coordinates will be updated to the target assembly)
GTF	GTF (Genome coordinates will be updated to the target assembly)
VCF	VCF (Genome coordinates and reference alleles will be updated to the target assembly)

#### SEVEN

### **USAGE**

#### Run CrossMap.py without any arguments will print help message:

```
# run CrossMap without argument
$ python CrossMap.py
Screen output:
Program: CrossMap (v0.1.1)
Description:
 CrossMap is a program for convenient conversion of genome coordinates
 and genomeannotation files between assemblies (eg. lift from human
 hg18 to hg19 or vice versa). It support file in BAM, SAM, BED, Wiggle,
 BigWig, GFF, GTF, VCF, etc.
Usage: CrossMap.py <command> [options]
 bam convert alignment file in BAM or SAM format.
 bed convert genome cooridnate or annotation file in BED or BED-like format.
            convert genome coordinate file in BigWig format.
 gff convert genome cooridnate or annotation file in GFF or GTF format.
 vcf convert genome coordinate file in VCF format.
 wig convert genome coordinate file in Wiggle, or bedGraph format.
```

#### Run CrossMap.py with command keyword will print help message for that command:

```
$ python CrossMap.py bed
```

#### Screen output:

```
Usage:
```

```
CrossMap.py bed input_chain_file input_bed_file [output_file]
```

#### Description:

```
"input_chain_file" and "input_bed_file" can be regular or compressed (*.gz, *.Z, *.z, *.bz, *.bz2, *.bzip2) file, local file or URL (http://, https://, ftp://) pointing to remote file. BED file must have at least 3 columns (chrom, start, end) and no more than 12 columns. If no "output_file" was specified, output will be directed to screen (console). BED format: http://genome.ucsc.edu/FAQ/FAQformat.html#format1
```

#### Example:

```
CrossMapy.py bed hg18ToHg19.over.chain.gz test.hg18.bed test.hg19.bed
# write output to "test.hg19.bed"
```

```
Example:
   CrossMapy.py bed hg18ToHg19.over.chain.gz test.hg18.bed
# write output to screen
```

#### 7.1 Convert BED format files

A BED (Browser Extensible Data) file is a tab-delimited text file describing genome regions or gene annotations. It is the standard file format used by UCSC. It consists of one line per feature, each containing 3-12 columns. CrossMap converts BED files with less than 12 columns to a different assembly by updating the chromosome and genome coordinates only; all other columns remain unchanged. Regions from old assembly mapping to multiple locations to the new assembly will be split. For 12-columns BED files, all columns will be updated accordingly except the 4th column (name of bed line), 5th column (score value) and 9th column (RGB value describing the display color). 12-column BED files usually define multiple blocks (eg. exon); if any of the exons fails to map to a new assembly, the whole BED line is skipped.

The input BED file can be plain text file, compressed file with extension of .gz, .Z, .z, .bz, .bz2 and .bzip2, or even a URL pointing to accessible remote files (http://, https:// and ftp://). Compressed remote files are not supported. The output is a BED format file with exact the same number of columns as the original one.

Standard BED format has 12 columns, but CrossMap also supports BED-like formats:

- BED3: The first 3 columns ("chrom", "start", "end") of BED format file.
- BED6: The first 6 columns ("chrom", "start", "end", "name", "score", "strand") of BED format file.
- Other: Format has at least 3 columns ("chrom", "start", "end") and no more than 12 columns. All other columns are arbitrary.

#### NOTE:

- 1. For BED-like formats mentioned above, CrossMap only updates "chrom (1st column)", "start (2nd column)", "end (3rd column)" and "strand" (if any). All other columns will keep AS-IS.
- 2. Lines starting with '#', 'browser', 'track' will be skipped.
- 3. Lines will less than 3 columns will be skipped.
- 4. 2nd-column and 3-column must be integer, otherwise skipped.
- 5. "+" strand is assumed if no strand information was found.
- 6. For standard BED format (12 columns). If any of the defined exon blocks cannot be uniquely mapped to target assembly, the whole entry will be skipped.
- 7. "input\_chain\_file" and "input\_bed\_file" can be regular or compressed (.gz, .Z, .z, .bz, .bz2, .bzip2) file, local file or URL (http://, https://, ftp://) pointing to remote file.
- 8. If output\_file was not specified, results will be printed to screen (console). In this case, the original bed entries (include items failed to convert) were also printed out.
- 9. If input region cannot be consecutively mapped target assembly, it will be split.

Example (run CrossMap with **no** *output\_file* specified):

```
$ python CrossMap.py bed hg18ToHg19.over.chain.gz test.hg18.bed3
```

Conversion results were printed to screen directly (column1-3 are hg18 based, column5-7 are hg19 based):

16 Chapter 7. Usage

chr1	142614848	142617697	->	chr1	143903503	143906352
chr1	142617697	142623312	->	chr1	143906355	143911970
chr1	142623313	142623350	->	chr1	143911971	143912008
chr1	142623351	142626523	->	chr1	143912009	143915181
chr1	142633862	142633883	->	chr1	143922520	143922541
chr1	142633884	142636152	->	chr1	143922542	143924810
chr1	142636152	142636326	->	chr1	143924813	143924987
chr1	142636339	142636391	->	chr1	143925000	143925052
chr1	142636392	142637362	->	chr1	143925052	143926022
chr1	142637373	142639738	->	chr1	143926033	143928398
chr1	142639739	142639760	->	chr1	143928399	143928420
chr1	142639761	142640145	->	chr1	143928421	143928805
chr1	142640153	142641149	->	chr1	143928813	143929809

#### Example (run CrossMap with output\_file (test.hg19.bed3) specified):

\$ python CrossMap.py bed hg18ToHg19.over.chain.gz test.hg18.bed3 test.hg19.bed3

```
$ cat test.hg19.bed3
     143903503
chr1
                       143906352
chr1
      143906355
                       143911970
chr1
      143911971
                       143912008
      143912009
chr1
                       143915181
chr1
      143922520
                       143922541
chr1
      143922542
                       143924810
      143924813
chr1
                       143924987
      143925000
chr1
                       143925052
chr1
       143925052
                       143926022
chr1
       143926033
                       143928398
chr1
       143928399
                       143928420
chr1
       143928421
                       143928805
       143928813
                       143929809
chr1
```

#### Example (one input region was split because it cannot be consecutively mapped target assembly):

\$ python CrossMap.py bed hg18ToHg19.over.chain.gz test.hg18.bed3

chr10	81346644	81349952	+	->	chr10	81356692	81360000	0	+
chr10	81349952	81364937	+	->	chr10	81360000	81374985	5	+
chr10	81364952	81365854	+	->	chr10	81375000	81375902	2	+
chr10	81365875	81369946	+	->	chr10	81375929	81380000	0	+
chr10	81369946	81370453	+	->	chr10	81380000	8138050	7	+
chr10	81370483	81371363	+	->	chr10	81380539	81381419	9	+
chr10	81371363	81371365	+	->	chr10	62961832	6296183	4	+
chr10	81371412	81371432	+	(split	.1:chr10:	:81371412:8137142	.2:+)	chr10	62961
chr10	81371412	81371432	+	(split	.2:chr10:	:81371422:8137143	.2:+)	chrX	632783

#### 7.2 Convert BAM/SAM format files

SAM (Sequence Alignment Map) format is a generic format for storing sequencing alignments, and BAM is binary and compressed version of SAM (Li et al., 2009). Most high-throughput sequencing (HTS) alignments were in SAM/BAM format and many HTS analysis tools work with SAM/BAM format. CrossMap updates chromosomes, genome coordinates, header sections, and all SAM flags accordingly. The program version (of CrossMap) is inserted into the header section, along with the names of the original BAM file and the chain file. For pair-end sequencing, insert size is also recalculated. The input BAM file should be sorted and indexed properly using samTools (Li et al., 2009). Output format is determined from the input format and BAM output will be sorted and indexed automatically.

#### Typing command without any arguments will print help message:

```
$ python CrossMap.py bam
```

#### Screen output:

#### Example (Convert BAM from hg19 to hg18):

```
$ python2.7 CrossMap.py bam hg19ToHg18.over.chain.gz test.hg19.bam test.hg18
@ 2013-11-15 14:08:01: Read hg19ToHg18.over.chain.gz ...
@ 2013-11-15 14:08:01: Liftover BAM file: test.hg19.bam ==> test.hg18.bam
@ 2013-11-15 14:08:17: Done!
@ 2013-11-15 14:08:17: Total entries: 164930
@ 2013-11-15 14:08:17: Failed to map: 5257
@ 2013-11-15 14:08:17: Sort "test.hg18.bam" ...
@ 2013-11-15 14:08:23: Index "test.hg18.sorted.bam" ...
```

#### # BAM/SAM header sections was updated:

\$ samtools view -H test.hg18.bam

```
$ samtools view -H test.hg19.bam
    SN:chr1 LN:249250621
@SQ
@SQ
      SN:chr2 LN:243199373
@SO
      SN:chr3 LN:198022430
@SO
      SN:chr4 LN:191154276
@SO
    SN:chr5 LN:180915260
@SO
    SN:chr6 LN:171115067
@SQ
    SN:chr7 LN:159138663
@SQ
    SN:chr8 LN:146364022
    SN:chr9 LN:141213431
@SO
@SQ
     SN:chr10
                    LN:135534747
@SO
     SN:chr11
                    LN:135006516
@SQ
     SN:chr12
                    LN:133851895
@SQ
      SN:chr13
                     LN:115169878
@SQ
      SN:chr14
                     LN:107349540
@SQ
      SN:chr15
                     LN:102531392
@SO
      SN:chr16
                     LN:90354753
@SO
      SN:chr17
                     LN:81195210
@SQ
      SN:chr18
                     LN:78077248
@SO
      SN:chr19
                     LN:59128983
@SQ
     SN:chr20
                     LN: 63025520
@SQ
    SN:chr21
                     LN:48129895
@SO
    SN:chr22
                    LN:51304566
    SN:chrX LN:155270560
@SQ
@SO
    SN:chrY LN:59373566
     SN:chrM LN:16571
0.50
      ID:Sample_618545BE
                             SM:Sample_618545BE LB:Sample_618545BE PL:Illumina
arg
@PG
     ID:bwa PN:bwa VN:0.6.2-r126
```

18 Chapter 7. Usage

```
@HD
      VN:1.0 SO:coordinate
@SQ
      SN:chr1 LN:247249719
@SQ
      SN:chr10
                    LN:135374737
@SQ
      SN:chr11
                     LN:134452384
@SQ
      SN:chr11_random LN:215294
@SQ
      SN:chr12
                    LN:132349534
@SQ
      SN:chr13
                     LN:114142980
      SN:chr13_random LN:186858
@SQ
      SN:chr14 LN:106368585
@SQ
      SN:chr15
@SQ
                    LN:100338915
@SQ
      SN:chr15_random LN:784346
@SQ
      SN:chr16
                    LN:88827254
@SQ
      SN:chr17
                    LN:78774742
@SO
     SN:chr17_random LN:2617613
@SQ
      SN:chr18 LN:76117153
@SQ
      SN:chr18_random LN:4262
@SQ
      SN:chr19 LN:63811651
@SQ
      SN:chr19_random LN:301858
@SQ
      SN:chrl_random LN:1663265
      SN:chr2 LN:242951149
@SQ
                LN:62435964
@SO
      SN:chr20
@SQ
      SN:chr21
                     LN:46944323
@SQ
    SN:chr21_random LN:1679693
@SQ
    SN:chr22 LN:49691432
@SQ
    SN:chr22_random LN:257318
    SN:chr3 LN:199501827
@SQ
@SO
      SN:chr3 random LN:749256
@SQ
      SN:chr4 LN:191273063
@SQ
      SN:chr4_random LN:842648
@SQ
      SN:chr5 LN:180857866
@SQ
      SN:chr6 LN:170899992
@SO
      SN:chr6_random LN:1875562
@SQ
      SN:chr7 LN:158821424
@SO
      SN:chr7_random LN:549659
@SO
      SN:chr8 LN:146274826
@SQ
      SN:chr8_random LN:943810
@SQ
      SN:chr9 LN:140273252
@SQ
      SN:chr9_random LN:1146434
@SO
      SN:chrM LN:16571
@SO
      SN:chrX LN:154913754
aso
      SN:chrX_random LN:1719168
@SO
      SN:chrY LN:57772954
@RG
     ID:Sample_618545BE
                             SM:Sample_618545BE LB:Sample_618545BE PL:Illumina
@PG
      PN:bwa ID:bwa VN:0.6.2-r126
@PG
      ID:CrossMap
                     VN:0.1.3
a CO
      Liftover from original BAM/SAM file: test.hg19.bam
@CO
      Liftover is based on the chain file: ../test/hg19ToHg18.over.chain.gz
```

- 1. Input is BAM or SAM format file. Output format depends on input format. (i.e BAM -> BAM, SAM -> SAM)
- 2. Alignments that are failed to convert will be saved in ".unmap.bam" or '.unmap.sam'.
- 3. Header section will be updated to target assembly.
- 4. Genome coordinates and all SAM flags in alignment section will be updated to target assembly.
- 5. Optional fields in alignment section will not be updated in current version (v0.1.3).

# 7.3 Convert Wiggle/BigWig format files

Wiggle (WIG) format is useful for displaying continuous data such as GC content and reads intensity of high-throughput sequencing data. BigWig is a self-indexed binary-format Wiggle file, and has the advantage of supporting random access. This means only regions that need to be displayed are retrieved by genome browser, and it dramatically reduces the time needed for data transferring (Kent et al., 2010). Input wiggle data can be in variableStep (for data with irregular intervals) or fixedStep (for data with regular intervals). Regardless of the input, the output will always in bedGraph format bedGraph format is similar to wiggle format and can be converted into BigWig format using UCSC wigToBigWig tool. We export files in bedGraph because it is usually much smaller than file in wiggle format, and more importantly, CrossMap internally transforms wiggle into bedGraph to increase running speed.

If an input file is in BigWig format, the output is BigWig format if UCSC's 'wigToBigWig' executable can be found; otherwise, the output file will be in bedGraph format.

Typing command without any arguments will print help message:

```
$ python2.7 CrossMap.py wig

Screen output:

Usage:
    CrossMap.py wig input_chain_file input_wig_file output_prefix

Description:
    "input_chain_file" can be regular or compressed (*.gz, *.Z, *.z, *.bz, *.bz2, *.bzip2) file, local file or URL (http://, https://, ftp://) pointing to remote file. Both "variableStep" and "fixedStep" wiggle lines are supported. Wiggle format: http://genome.ucsc.edu/goldenPath/help/wiggle.html

Example:
    CrossMapy.py wig hg18ToHg19.over.chain.gz test.hg18.wig test.hg19
```

NOTE:

1. To improve performance, this script calls GNU "sort" command internally. If "sort" command does not exist, CrossMap will exit.

Typing command without any arguments will print help message:

```
$ python2.7 CrossMap.py bigwig

Screen output:

Usage:
    CrossMap.py bigwig input_chain_file input__bigwig_file output_prefix

Description:
    "input_chain_file" can be regular or compressed (*.gz, *.Z, *.z, *.bz, *.bz2, *.bzip2) file, local file or URL (http://, https://, ftp://) pointing to remote file. Bigwig format: http://genome.ucsc.edu/goldenPath/help/bigWig.html

Example:
    CrossMapy.py bigwig hg18ToHg19.over.chain.gz test.hg18.bw test.hg19

Example (Convert BigWig file from hg18 to hg19):

$ python CrossMap.py bigwig hg19ToHg18.over.chain.gz test.hg19.bw test.hg18

@ 2013-11-17 22:12:42: Read chain_file: ../data/hg19ToHg18.over.chain.gz

@ 2013-11-17 22:12:44: Liftover bigwig file: test.hg19.bw ==> test.hg18.bgr

@ 2013-11-17 22:15:38: Merging overlapped entries in bedGraph file ...
```

20 Chapter 7. Usage

```
@ 2013-11-17 22:15:38: Sorting bedGraph file:test.hg18.bgr
@ 2013-11-17 22:15:39: Convert wiggle to bigwig ...
```

- 1. To improve performance, this script calls GNU "sort" command internally. If "sort" command does not exist, CrossMap will exit.
- 2. Output files: output\_prefix.bw, output\_prefix.bgr, output\_prefix.sorted.bgr

#### 7.4 Convert GFF/GTF format files

GFF (General Feature Format) is another plain text file used to describe gene structure. GTF (Gene Transfer Format) is a refined version of GTF. The first eight fields are the same as GFF. Plain text, compressed plain text, and URLs pointing to remote files are all supported. Only chromosome and genome coordinates are updated. The format of output is determined from the input.

```
Typing command without any arguments will print help message:

$ python2.7 CrossMap.py gff

Screen output:

Usage:
    CrossMap.py gff input_chain_file input_gff_file output_file
```

```
Description:

"input_chain_file" can be regular or compressed (*.gz, *.Z, *.z, *.bz, *.bz2, *.bzip2) file, local file or URL (http://, https://, ftp://) pointing to remote file. input file must be in GFF or GTF format. GFF format: http://genome.ucsc.edu/FAQ/FAQformat.html#format3 GTF format: http://genome.ucsc.edu/FAQ/FAQformat.html#format4
```

#### Example:

rample.

CrossMap.py gff hg19ToHg18.over.chain.gz test.hg19.gtf test.hg18.gtf #write output to test.hg18.gr

#### Example:

CrossMap.py gff hg19ToHg18.over.chain.gz test.hg19.gtf # write output to screen

\$ python CrossMap.py gff hg19ToHg18.over.chain.gz test.hg19.gtf test.hg18.gtf

#### Example (Convert GTF file from hg19 to hg18):

```
@ 2013-11-17 20:44:47: Read chain_file: ../data/hg19ToHg18.over.chain.gz
$ head test.hg19.gtf
chr1 hg19_refGene
                     CDS
                             48267145
                                             48267291
                                                            0.000000
                                                                                   ()
                                                                                           gene_:
chr1 hg19_refGene
                             66081691
                                             66081907
                                                            0.000000
                     exon
                                                                                           gene_:
chr1 hg19_refGene CDS
                             145334684
                                             145334792
                                                            0.000000
                                                                                   2.
                                                                                           gene_:
chr1 hg19_refGene exon
                             172017752
                                            172017890
                                                            0.000000
                                                                                           gene_:
chrl hg19_refGene CDS
                             206589249
                                             206589333
                                                            0.000000
                                                                                   2.
                                                                                           gene :
```

```
chr1 hg19_refGene exon 210573812
                                         210574006
                                                       0.000000
                                                                                    gene_:
chrl hg19_refGene CDS
                          235850249
                                         235850347
                                                       0.000000
                                                                                    gene_:
chr1 hg19_refGene CDS
                          235880012
                                        235880078
                                                       0.000000
                                                                            1.
                                                                                    gene
chrl hg19_refGene exon 3417741 3417872 0.000000
                                                                     gene_id "NM_001409";
    hg19_refGene
                                         10190871
                                                       0.000000
chr1
                   exon
                          10190773
                                                                                    gene_:
$
```

\$ head	d test.hg18.gtf							
chr1	hg19_refGene	CDS	48039732	48039878	0.000000	_	0	gene_:

```
65854279
                                              65854495
                                                               0.000000
chr1
     hg19_refGene
                                                                                              gene_:
                      exon
                                                              0.000000
      hg19_refGene
                      CDS
                              144046041
                                              144046149
                                                                                      2
chr1
                                                                                              gene_:
      hg19_refGene
                              170284375
                                              170284513
                                                               0.000000
chr1
                      exon
                                                                                              gene_:
chr1
      hg19_refGene
                      CDS
                              204655872
                                              204655956
                                                               0.000000
                                                                                      2
                                                                                              gene_:
      hg19_refGene
                      exon
                              208640435
                                              208640629
                                                               0.000000
                                                                                              gene_:
                                                              0.000000
chr1
      hg19_refGene
                      CDS
                              233916872
                                              233916970
                                                                                              gene_:
     hg19_refGene
                                                                                              gene_:
chr1
                      CDS
                              233946635
                                              233946701
                                                              0.000000
                                                                                      1
                              3407601 3407732 0.000000
                                                                              gene_id "NM_001409";
chr1 hg19_refGene
                      exon
                                              10113458
                                                              0.000000
chr1 hg19_refGene
                              10113360
                      exon
                                                                                              gene_:
```

- 1. Each feature (exon, intron, UTR, etc) is processed separately and independently, and we do NOT check if features originally belonging to the same gene were converted into the same gene.
- 2. If user want to liftover gene annotation files, use BED12 format.
- 3. If no output file was specified, output will be printed to screen (console). In this case, items failed to convert are also printed out.

#### 7.5 Convert VCF format files

\$ python2.7 CrossMap.py qff

VCF (variant call format) is a flexible and extendable line-oriented text format developed by the 1000 Genome Project. It is useful for representing single nucleotide variants, indels, copy number variants, and structural variants. Chromosomes, coordinates, and reference alleles are updated to a new assembly, and all the other fields are not changed.

Typing command without any arguments will print help message:

```
Screen output:

usage:
    CrossMap.py vcf input_chain_file input_VCF_file ref_genome_file output_file

Description:
    "input_chain_file" and "input_VCF_file" can be regular or compressed (*.gz, *.Z, *.z, *.bz, *.bz2, *.bzip2) file, local file or URL (http://, https://, ftp://) pointing to remote file. "ref_genome_file" is genome sequence file of 'target assembly' in FASTA foramt.

Example:
    CrossMap.py vcf hg19ToHg18.over.chain.gz test.hg19.vcf hg18.fa test.hg18.vcf
```

#### Example (Convert VCF file from hg19 to hg18):

```
$ python CrossMap.py vcf hg19ToHg18.over.chain.gz test.hg19.vcf ../database/genome/hg18.fa test.hg19
@ 2013-11-17 20:53:37: Read chain_file: ../data/hg19ToHg18.over.chain.gz
@ 2013-11-17 20:53:39: Load reference genome: ../../database/genome/hg19.fa
cannot fetch sequence from ../../database/genome/hg19.fa for chr21_random:1363681-1363682
@ 2013-11-17 20:57:00: Total entries: 998
@ 2013-11-17 20:57:00: Failed to map: 2
$ grep -v '#' test.hg19.vcf |head -10
chr1
     10933566
                              С
                                      G
                                                       PASS
                                                               ADP=13; WT=0; HET=0; HOM=1; NC=0
                                                                                               GT:GQ
     11187893
                              Τ
                                      С
                                                       PASS
                                                               ADP=224; WT=0; HET=0; HOM=1; NC=0
                                                                                               GT:GO
chr1
chr1 11205058
                              С
                                      Т
                                                       PASS
                                                               ADP=625; WT=0; HET=0; HOM=1; NC=0
                                                                                               GT:GQ
                                       G
                                                       PASS
                                                               ADP=52; WT=0; HET=0; HOM=1; NC=0
chr1 11292753
                              Α
                                                                                               GT:GQ
```

22 Chapter 7. Usage

chr1	11318763		С	G		str10	ADP=88; WT=0; HET=0; HOM=1; NC=0	GT:GQ
chr1	11319587		A	G		PASS	ADP=70; WT=0; HET=0; HOM=1; NC=0	GT:GQ
chr1	16202995		С	T		PASS	ADP=463; WT=0; HET=1; HOM=0; NC=0	GT:GQ
chr1	27088546	•	А	T		PASS	ADP=124; WT=0; HET=1; HOM=0; NC=0	GT:GQ:
chr1	27101390		T	С		str10	ADP=267; WT=0; HET=1; HOM=0; NC=0	GT:GQ
chr1	34007097	•	Т	С	•	PASS	ADP=10;WT=0;HET=1;HOM=0;NC=0	GT:GQ
\$ grep	-v '#' test.hg1	.8.vcf	head -10	)				
1	10856153		С	G		PASS	ADP=13;WT=0;HET=0;HOM=1;NC=0	GT:GQ
1	11110480		T	С		PASS	ADP=224; WT=0; HET=0; HOM=1; NC=0	GT:GQ
1	11127645		С	T	•	PASS	ADP=625; WT=0; HET=0; HOM=1; NC=0	GT:GQ
1	11215340		A	G		PASS	ADP=52; WT=0; HET=0; HOM=1; NC=0	GT:GQ
1	11241350		С	G		str10	ADP=88; WT=0; HET=0; HOM=1; NC=0	GT:GQ
1	11242174		A	G		PASS	ADP=70; WT=0; HET=0; HOM=1; NC=0	GT:GQ
1	16075582		С	T	•	PASS	ADP=463;WT=0;HET=1;HOM=0;NC=0	GT:GQ
1	26961133		A	T	•	PASS	ADP=124; WT=0; HET=1; HOM=0; NC=0	GT:GQ
1	26973977		T	С		str10	ADP=267; WT=0; HET=1; HOM=0; NC=0	GT:GQ
1	33779684	•	T	С	•	PASS	ADP=10;WT=0;HET=1;HOM=0;NC=0	GT:GQ
\$ grep	-v '#' test.hg1	.8.vcf.u	ınmap	#coordi	nates ar	re still	based on hg19	
chr14	20084444		G	С		PASS	ADP=253; WT=0; HET=1; HOM=0; NC=0	GT:GQ
chr14	20086290		T	С		PASS	ADP=441; WT=0; HET=1; HOM=0; NC=0	GT:GQ:

- 1. Genome coordinates and reference allele will be updated to target assembly.
- 2. Reference genome is genome sequence of target assembly.
- 3. Output files: output\_file and output\_file.unmap.

24 Chapter 7. Usage


# COMPARE TO UCSC LIFTOVER TOOL

To access the accuracy of CrossMap, we randomly generated 10,000 genome intervals (download from here) with the fixed interval size of 200 bp from hg19. Then we converted them into hg18 using CrossMap and UCSC liftover tool with default configurations. We compare CrossMap to UCSC liftover tool because it is the most widely used tool to convert genome coordinates.

CrossMap failed to convert 613 intervals, and UCSC liftover tool failed to convert 614 intervals. All failed intervals were exactly the same except one region (chr2 90542908 90543108). UCSC failed to convert it because this region needs to be split 2 times:


Original (hg19)	Split (hg19)	Target (hg18)
chr2 90542908 90543108 -	chr2 90542908 90542933 -	chr2 89906445 89906470 -
chr2 90542908 90543108 -	chr2 90542933 90543001 -	chr2 87414583 87414651 -
chr2 90542908 90543108 -	chr2 90543010 90543108 -	chr2 87414276 87414374 -

For genome intervals that were successfully converted to hg18, the start and end coordinates were exactly the same between UCSC conversion and CrossMap conversion.



Start coordinates (UCSC liftover)

Start coordinates



End coordinates

25

#### **CHAPTER**

# **NINE**

#### **LICENSE**

#### CrossMap is distributed under GNU General Public License

This program is free software; you can redistribute it and/or modify it under the terms of the GNU General Public License as published by the Free Software Foundation; either version 2 of the License, or (at your option) any later version. This program is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU General Public License for more details. You should have received a copy of the GNU General Public License along with this program; if not, write to the Free Software Foundation, Inc., 51 Franklin Street, Fifth Floor, Boston, MA 02110-1301 USA

28 Chapter 9. LICENSE

# **CONTACT**

• Wang.Liguo AT mayo.edu





Making Cancer History®

30

